

20L AQUAPure[®]

Water for Injection

Sterile. Reliable. Ready for Integration.

Engineered Tray Designed for Cleanroom/Lab Efficiency & Workflow Versatility

Aquability offers complete Single Use System integration, providing a streamlined, efficient, and contamination-resistant solution for your process. Our ready-to-use **20L Water for Injection (WFI)** tray system reduces contamination risks, enhances operational speed, and supports cost-effective, regulatory-compliant production. Backed by a vertically integrated supply chain and a dedicated Quality Assurance Team, Aquability is your trusted partner in advancing Single Use System integration with confidence.



i For Additional Solutions
(**AQUA**Lutions), Custom Packaging,
and Container Sizes reach out to us
at Customercare@aquability.com

Key Features and Benefits



Unmatched Single-Use System Expertise

- 🔹 **Extensive Experience** in Single-Use Manufacturing and package engineering
- 🔹 **Complimentary Support** for seamless implementation
- 🔹 **Expert Package Engineering** enhances both bag and integrated product tray for cost-effective, high-efficiency production
- 🔹 **Fully Integrated In-House** bag manufacturing and filling for total quality control



Packaging That Performs

- 🔹 **Integrated Tray and Bag System** enhances contamination control and simplifies handling
- 🔹 **Poly Overwrap** eliminates cardboard use, minimizing particulate risk
- 🔹 **Invertible Lid Design** enables complete fluid evacuation when rotated 180°
- 🔹 **Durable Packaging** ensures leak-free delivery every time



High Performance Tubing Built In

- 🔹 **3/8" Wall Tubing** housed in an internal chase for cleanliness and tangle-free use
- 🔹 **Optimized Tubing Format** includes thick wall tubing ideal for peristaltic pump dispensing and additional welding and sealing
- 🔹 **Includes Sterile Male and Female Quick Connectors**, securely covered to ensure cleanliness



Environmentally Conscious Design

- 🔹 **Durable Trays** are designed for secondary use in operational workflows and may even be repurposed for final product packaging, extending their life cycle beyond the initial application
- 🔹 **Constructed from Fully Recyclable PETG** to support sustainable, eco-friendly practices

Contact us to learn more

<https://www.aquabiliti.com/contact-us>

or contact your Aquabiliti Account Manager

AQUABILITI[™]
IMPACTING LIVES



For Tray Evacuation Demonstration See [Here](#)



Aquabiliti Customer Service (615) 833-2633
www.aquabiliti.com

© Copyright 2024 Aquabiliti Life Sciences, LLC. All rights reserved.

YRK-MKT-MKT-0023 Rev. 1