

AQUAPure HYPER

Water for Injection (WFI) USP/EP Grade

Aquability's **AQUAPure HYPER Water for Injection (WFI) USP/EP Grade** enables our clients to refine and improve workflows to enhance productivity – from R&D, pilot to commercial production – by providing WFI quality water that complies with compendial (USP/EP) requirements.



Multiple Packaging Sizes Available



Readily Available



Compliant With Compendial (USP/EP) Specifications and Requirements



Custom Manufacturing Available



Leverage Our Continuous Manufacturing Capabilities to Your Advantage

AQUAPure HYPER	PACKAGED PRODUCT ANALYSIS	METHOD OF REFERENCE	SPECIFICATION	
	Appearance	Visual	Clear & Colorless Liquid	
	Total Organic Carbon	USP<643>	< 8000 ppb	
	Conductivity	USP<645>	≤ 5 μS/cm	
	Endotoxin	USP<85>	≤ 0.05 EU/mL ¹	
	Sterility	USP<71>	No Growth	
	Particulate Matter	USP<788>	≥ 25 μm in size	≤ 3 particles/mL
			≥ 10 μm in size	≤ 25 particles/mL
	Oxidizable Substances	EP Sterilized WFI	Solution remains faintly pink	
	Residue on Evaporation	EP Sterilized WFI	≤ 3 mg (0.003%)	
	Nitrates	EP Sterilized WFI	≤ 0.2 ppm	


¹ The USP Specification for Compendial Water requires endotoxin levels to be less than 0.25EU/mL. Aquability applies a tightened specification, limiting endotoxin to be less than or equal to 0.05 EU/mL.

WFI Generation Online Analysis	AQUAPure	WFI-USP
Conductivity	≤ 0.06 μS/cm	≤ 1.3 μS/cm
TOC	≤ 20 ppb	< 500 ppb

Aquability's AQUAPure HYPER Water for Injection (WFI) USP/EP Grade is available in multiple packaging options.



1,000L



200L



100L



50L



20L



10L

Aquability's AQUAPure HYPER Water for Injection (WFI) USP/EP Grade

- Meets the USP/EP monograph for Water for Injection packaged in bulk for commercial use
- Catalog made-to-stock products*
- Custom sizes/packaging available on-request

** Availability dependent on inventory at time of order. Prioritized made-to-stock (MTS) items include 20L and 200L package sizes. Other package sizes may not be consistently MTS and may be made-to-order (MTO) at time of customer order.*

Contact us to learn more

<https://www.aquability.com/contact-us>

or contact your Aquability Account Manager

AQUABILITY™
IMPACTING LIVES

Aquability Customer Service (615) 833-2633
www.aquability.com

© Copyright 2024 Aquability Life Sciences, LLC. All rights reserved.

