

# AQUAPure ULTRA

## Water for Injection (WFI) USP Grade

Aquabilti's **AQUAPure ULTRA Water for Injection (WFI) USP Grade** enables our clients to refine and improve workflows to enhance productivity – from R&D, pilot to commercial production – by providing WFI quality water that complies with compendial (USP) requirements.



Multiple Packaging  
Sizes Available



Readily  
Available



Compliant With  
Compendial (USP)  
Specifications and  
Requirements



Custom  
Manufacturing  
Available



Leverage Our  
Continuous  
Manufacturing  
Capabilities to  
Your Advantage

AQUAPure ULTRA	PACKAGED PRODUCT ANALYSIS	METHOD OF REFERENCE	SPECIFICATION	
	Appearance	Visual	Clear & Colorless Liquid	
	Total Organic Carbon	USP<643>	< 8000 ppb	
	Conductivity	USP<645>	≤ 5 μS/cm	
	Particulate Matter	USP<788>	≥ 25 μm in size	≤ 3 particles/mL
			≥ 10 μm in size	≤ 25 particles/mL
	Endotoxin	USP<85>	≤ 0.05 EU/mL <sup>1</sup>	
	Sterility	USP<71>	No Growth	


<sup>1</sup> The USP specification for Compendial Water requires endotoxin levels to be less than 0.25EU/mL. Aquabilti applies a tightened specification, limiting endotoxin to less than or equal to 0.05 EU/mL.

WFI Generation Online Analysis	AQUAPure	WFI-USP
Conductivity	≤ 0.06 μS/cm	≤ 1.3 μS/cm
TOC	≤ 20 ppb	< 500 ppb

# Aquabiliti's AQUAPure ULTRA Water for Injection (WFI) USP Grade is available in multiple packaging options.




1,000L



200L



100L



50L



20L



10L

## Aquabiliti's AQUAPure ULTRA Water for Injection (WFI) USP Grade

- Meets the USP monograph for Water for Injection packaged in bulk for commercial use
- Catalog made-to-stock products\*
- Custom sizes/packaging available on-request

*\* Availability dependent on inventory at time of order. Prioritized made-to-stock (MTS) items include 20L and 200L package sizes. Other package sizes may not be consistently MTS and may be made-to-order (MTO) at time of customer order.*

**Contact us to learn more**

<https://www.aquabiliti.com/contact-us>

or contact your Aquabiliti Account Manager



**AQUABILITI™**  
IMPACTING LIVES

**Aquabiliti Customer Service (615) 833-2633**  
[www.aquabiliti.com](http://www.aquabiliti.com)

© Copyright 2024 Aquabiliti Life Sciences, LLC. All rights reserved.