

# AQUAPure WIR

## Water for Irrigation Quality Water USP Grade

Aquabilti's **AQUAPure WIR Water for Irrigation Quality Water USP Grade** is a sterile non-pyrogenic water produced by our proprietary water purification system. **AQUAPure WIR** is ideal for use in the manufacture and testing of medical devices, and applications where Water for Injection may not be appropriate.



Multiple Packaging Sizes Available



Readily Available



Compliant With Compendial (USP) Specifications and Requirements



Custom Manufacturing Available



Leverage Our Continuous Manufacturing Capabilities to Your Advantage

AQUAPure WIR	PACKAGED PRODUCT ANALYSIS	METHOD OF REFERENCE	SPECIFICATION
	Appearance	Visual	Clear & Colorless Liquid
	Total Organic Carbon	USP<643>	< 8000 ppb
	Conductivity	USP<645>	≤ 5 μS/cm
	Endotoxin	USP<85>	≤ 0.05 EU/mL <sup>1</sup>
	Sterility	USP<71>	No Growth

<sup>1</sup> The USP specification for Compendial Water requires endotoxin levels to be less than 0.25 EU/mL. Aquabilti applies a tightened specification, limiting endotoxin to less than or equal to 0.05 EU/mL.

# Aquability's AQUAPure WIR Water for Irrigation Quality Water USP Grade

is available in multiple  
packaging options.



1,000L



200L



100L



50L



20L



10L

## Aquability's AQUAPure WIR Water for Irrigation Quality Water USP Grade

- Meets the USP monograph for Water for Irrigation Quality Water packaged in bulk for commercial use
- Catalog made-to-stock products\*
- Custom sizes/packaging available on-request

*\* Availability dependent on inventory at time of order. Prioritized made-to-stock (MTS) items include 20L and 200L package sizes. Other package sizes may not be consistently MTS and may be made-to-order (MTO) at time of customer order.*

**Contact us to learn more**

<https://www.aquability.com/contact-us>

or contact your Aquability Account Manager



**AQUABILITY™**  
IMPACTING LIVES

Aquability Customer Service (615) 833-2633  
[www.aquability.com](http://www.aquability.com)

© Copyright 2024 Aquability Life Sciences, LLC. All rights reserved.